

AeroDawg 

Canine Aerosol Chamber




SMALL

MEDIUM

LARGE

Mask Sizer

Each device is shipped with 2 different sized masks.



AeroDawg* Chamber SMALL,
includes small and medium masks

AeroDawg* Chamber LARGE,
includes medium and large masks

How to Use

The mask selected should cover the nose and mouth. However, if not possible apply mask to the dog's nostrils while gently holding the muzzle closed.



Trudell Medical International

www.aerodawg.com



201092-001 Rev.A. * trade-marks and registered trade-marks of Trudell Medical International.
© TMI 2015. All rights reserved.

AeroKat
Feline Aerosol Chamber

MEDIUM

SMALL

The exclusive **Flow-Vu**^{*} Inhalation Indicator moves toward the cat or dog as they inhale and will only move if there is a good seal.



Mask Sizer

Each **AeroKat**^{*} Chamber is shipped with a small and medium mask. The mask selected should cover the cat's nose and mouth but not the eyes.

Helping animals around the world breathe easier

www.aerokat.com

